

# **OxiTop BOD Measurement Instrumentation**

Pressure sensors for BOD testing and measurements.



## **Specifications**

Simple operation

measurements

easy operation

Non-Toxic

Measuring principle: Respirometric (manometric)
Display: 2-digit 7-segment LED, 10 mm
Measured parameter: Digits

Measuring range: 0 to 40 digits (+10 digits for overranging)

Stores 100 user-defined routine

User-friendly on screen guidance for

Accuracy: <u>+</u> 0.03 psi (<u>+</u> 2 hPa) or <u>+</u> 1 digit

Admissible measuring temperature: 68°F ± 1.8°F (20°C ± 1°C)
Admissible sample temperature (when being filled): 59-68°F (15-20°C)
Power supply: Lithium batteries (260 mAh), 2 x type CR 2430 (3V)

Power consumption: Max. 25 mA (during measurement)

Battery run time: ~ 3 years with normal use

Protection class: 3 IEC 1010 Protection system: IP 54 IEC 529

EMC

Emissions: EN 50081-1, FCC Class A

Immunity: EN 50082-2, NAMUR recommendation

Climate class: 2, VDI/VDE 3540

Ambient temperature:

Storage: -13°F to +149°F (-25°C to +65°C) Operation: 68°F ± 1.8 F (20°C ± 1 C)

Test mark: CE

Dimensions: H: 2.7in, dia. 2.8in (69 mm, dia. 70 mm)

Weight: Approx. 3 oz. (85 g)

## **Options and Accessories**

#### 208214 OxiTop® IS 6

BOD Measurement Instrumentation

BOD measuring instrument for self-check with

BOD measuring instrument for self-check with 6 measuring units, 115 VAC, 50/60 Hz power plug, and accessories.

#### 208215 OxiTop® IS 12

BOD Measurement Instrumentation BOD measuring instrument for self-check with 12 measuring units, 115 VAC, 50/60 Hz power plug, and accessories.

#### 208433 OxiTop® Thermostat Box

Thermostat Box

OxiTop® BOD Instrumentation thermostat box, for up to 12 OxiTop® units and stirrer base, 115 VAC, 50 Hz.

**Regional Distributor** 

### **Description**

#### OxiTop IS 6, IS 12

Measurement using OxiTop® Biological Oxygen Demand (BOD) Instrumentation is based on pressure measurement in a closed system: microorganisms in the sample consume the oxygen and form CO<sup>2</sup>. This is absorbed by NaOH, creating a vacuum, which can be read directly as a measured value in mg/l BOD.

The sample volume being tested regulates the amount of oxygen available for a complete BOD measurement. BOD measurement ranges of up to 4,000 mg/l can be measured using different sample volumes.

The OxiTop® BOD instrumentation heads (green and yellow for differentiation of inflow/outflow) have an AutoTemp function: if the sample temperature is too cold, the start of BOD measurement is automatically delayed (by at least 1 hour) until a constant temperature has been reached.

Apart from the automatic storage of 5 measured values (1 value per day), further measured BOD values can be read at all times during or after the period of 5 days, which permits the tracking of check values or measurements over longer periods. The OxiTop BOD Instrumentation measurement system also offers simple operation, improved controllability, and non-toxicity.



Contact

needs:

**Global Water** 

instrumentation

for all your

Water Level



